
# Vibe

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
 Mein Held  
[doccheck.com](https://doccheck.com)

# Vibe

## User Manual

### Intended use

Hi there, I'm Vibe! I'm intended solely for use by medical professionals to assist in medical diagnostics.

I'm a pallometer and serve to assess vibration sensitivity, allowing you to detect impairments of the peripheral nervous system.

### Improper use

If you use me in a way that differs from the instructions in this manual, it is considered improper use.

Please note that I am not waterproof, so don't test my swimming skills.

The manufacturer accepts no liability for any damage resulting from improper use. In such cases, the risk lies entirely with you, the user. So please be careful!

### Accessories

- 1 USB-C cable
- 1 USB-C cover

### Environmental conditions

I am designed for use in professional healthcare facilities.

I am not suitable for use in highly electromagnetic environments, particularly in areas with strong electromagnetic fields (e.g. MRI rooms) or explosive zones.

### Operation

#### Check and functional test

Please carry out a visual inspection for damage and wear before each use. Please also check if I function correctly. Do not use me if you detect any issues!

Before use, check that my sleeve **4** is fully attached, otherwise vibration transmission cannot be guaranteed.

You can recognize this by the fact that it ends at my clip.

#### Charging the battery and usage time

I'm equipped with a rechargeable lithium battery. To charge me, connect the supplied USB-C cable **3** to a USB power adapter.

Only use CE-certified USB power adapters from reputable manufacturers that meet the safety requirements of regulation IEC 60601-1 (Medical Electrical Equipment).

I have a power consumption of 160–320 mA. Please do not use any adapters or USB-C cables that are visibly damaged.

My battery should be charged outside the patient environment (at least 1.5 m distance). Please fully charge me before first use using the supplied USB-C cable.

Once charging is complete, my blue indicator light will turn off.

I can be used continuously for up to 2 hours before needing to be charged again.

If my charge level is too low, my button **1** will flash red for 10 seconds. This is your sign to stop using me and charge me up again. Only this way can the battery voltage remain high enough to properly power the motor and achieve the correct vibration intensity.

I also have a sleep mode. If I'm turned on but not in use, I'll automatically enter sleep mode after 10 minutes. In this mode, my LED and motor are off, and my CPU will enter energy-saving mode. To wake me up, just press the button **1** or toggle the switch **2** off and on again.

#### Usage


Turn me on using the switch **2** on the side. I vibrate briefly and flash once quickly to show you that I am ready.

Now place my front end, including my sleeve **4**, on

prominent bony points like the kneecap or radial head.


Press my button **1** to start vibration on level 1. You'll know it's level 1 because my button will light up red.

This corresponds to scale value 2/8 on a tuning fork.


Repeat the test at different vibration levels. Proceed as follows:

Press the button again – my vibration will reduce, and the LED will change to yellow. This corresponds to 4/8 on the tuning fork scale.

Press a third time – the lowest vibration level will activate, the LED will turn green. This corresponds to scale value of 6/8.

Press a fourth time – vibration will stop; the LED will turn off.

Press again – vibration will start again at level 1 (2/8).

Once the test is complete, please turn me off using the switch **2**.

Please make sure to test several points and compare sides (left and right). If the patient doesn't feel any vibration, this may indicate a neurological impairment and further examinations should be considered.

Note: I do not use the colors red, yellow, and green to indicate system status or alerts. My color logic is used to help you clinically interpret the patient's vibration perception. The levels follow a traffic light system for intuitive assessment.

## Materials

I owe my elegant shape to my anodized aluminum housing. Its two parts are held together with stainless steel screws. Thanks to my stainless steel clip, you can attach me to your coat pocket.

My sleeve **4** is made of biocompatible silicone.

On the inside, I contain a PCB and an ICR 10440 lithium battery. I also have a motor built in to generate vibrations.

All my materials that may come into direct contact with skin are considered skin-compatible and very rarely cause allergic reactions. If you suspect intolerance or allergy to any of the materials used, I should be used cautiously or not at all.

## Replacing the battery


If my battery stops working, you can easily replace it yourself.

Before opening my housing, please make sure that I am switched off and that my motor is not moving!


Place me on a solid surface. Use a suitable Allen key to loosen the first screw. Now you can remove the clip and take out the second screw.


Take off the upper part of my housing.

Leave my circuit board, as well as the switch and the motor, in the housing. If the button or the indicator light move and bother you during the replacement, you can briefly put them aside.


Changing my battery is very easy. First, remember how the battery and its cables are positioned in the housing. Then remove the battery and pull the plug out of the connector. This is easiest if you hold the circuit board close to the connector with one hand and gently pull on the battery cable with the other hand.

Now insert a new battery using the same principle. The connection only fits one way, so you can't go wrong when plugging in the new battery.

In the end, place the battery back into my housing as before. Make sure the cables are laid around the screw fixtures as before, otherwise you won't be able to close my housing properly.

If you previously removed the button or the indicator light, put them back in, close my housing, and screw it back together.

Done. I am ready for use again.

## Cleaning, disinfection and storage

### Manual cleaning

Please clean me regularly. Doing so is very simple: Wipe me down thoroughly with a damp cloth. My manufacturer recommends wipes with 70 % isopropyl alcohol or soapy water.

If you use other cleaning/disinfection agents, test their compatibility with my materials first. Otherwise, my sleek appearance might suffer. Frequent cleaning doesn't bother me. I like being clean.


### Sterilization

Do NOT sterilize me. If you do, I will break and probably come to a sad end.


## Storage

I'm not picky. It is best to store me in a dry place, inside my packaging, and at room temperature. I have nothing against cool rooms, either.

My sleeve **4** has a usage position and a transport position. In the transport position, it sits over the button **1** to prevent accidental activation

Usage position

Transport position

## Maintenance and warranty

If my battery is defective, you can replace it. The warranty becomes void if replacement parts not purchased from DocCheck Shop are used.

As long as you use me correctly and follow this manual, the warranty covers manufacturing and material defects for two years starting on the date of purchase. Obvious misuse and accidental damage are excluded. Cosmetic defects such as surface scratches are also excluded from the warranty.

## Disposal

Please note that this product must be disposed of as separate electrical and electronic waste. If you need to dispose of the battery, take it to the appropriate collection site. Please consult your local disposal regulations for specifics.

## Safety notice

I am not intended for invasive use.

Due to small parts (USB-C cover), do not give me to children under 3 years!

Check that I function properly before each use. If damaged, stop using me immediately. If I no longer function, only allow repairs by qualified personnel.

Also, do not use me in fire- or explosion-hazardous environments or near strong magnetic fields.

### Note to users and/or patients

All serious incidents associated with me must be reported to DocCheck Steam GmbH at [qm@docchecksteam.com](mailto:qm@docchecksteam.com) and to the competent authority of the Member State.

## Electromagnetic safety

### Exclusions

I'm not designed for operation in special environments as defined in IEC standard 60601-1-2, including:

- Electrosurgery environments with high electromagnetic disturbances.
- Areas with strong electromagnetic fields, such as those generated by radiotherapy systems.

### Precautions to avoid electromagnetic interference

To avoid adverse events caused by electromagnetic interference, please observe the following instructions:

### Use of accessories

Only use accessories and cables that are specifically approved by my manufacturer for this device. The use of unauthorized accessories may increase electromagnetic emissions or reduce immunity to interference, which can lead to malfunctions.

### Distance from high-frequency devices

Portable devices that use high-frequency communication (e.g. Wi-Fi peripherals, antenna cables, or external antennas), including cables specified by my manufacturer, should be kept at least 30 cm (12 inches) away from me.

### Use near other devices

Avoid using me in close proximity to or stacked with other equipment. If it can't be avoided, please monitor how the devices behave when used together to ensure that both function properly.

### Protection against electrostatic discharge (ESD)

Electrostatic discharges can impair my functionality. Therefore, I should be used in environments that minimize electrostatic discharge (ESD), such as those with suitable flooring or controlled humidity.

Immunity testing	Basic EMC standard	EN/IEC 60601 Test conditions	Compliance
Magnetic fields at utility frequency (50/60 Hz)	EN/IEC 61000-4-8	30 A/m 50 and 60 Hz	30 A/m 50 and 60 Hz
Electrostatic discharge	IEC 61000-4-2	Contact Discharge: ±8 kV  Air Discharge: ±2 kV, ±4 kV, ±8 kV, ±15 kV	Contact Discharge: ±8 kV  Air Discharge: ±2 kV, ±4 kV, ±8 kV, ±15 kV
High-frequency electromagnetic fields	IEC 61000-4-3	3 V/m 80 MHz–2.7 GHz 80% AM at 1 kHz	3 V/m 80 MHz–2.7 GHz

### Measures in case of electromagnetic interference

If I cause or am affected by interference, the following actions can help ensure my basic safety:

- Reposition me or use me in a different location.
- Increase the distance between me and other electronic devices.
- Connect me to a different power source when charging.

### Essential performance

I do not possess essential performance characteristics that would pose an unacceptable risk to patients, users, or third parties in case of loss or impairment.

### Compliance with emission and immunity standards

I meet the requirements of IEC standard 60601-1-2 regarding:

- Emissions: CISPR 11 Class B, Group 1
- Immunity: see the following table:

### Manufacturer


DocCheck Steam GmbH  
Vogelsanger Str. 66, 50823 Cologne, Germany  
[doccheckshop.eu/more/doccheck-things/](https://doccheckshop.eu/more/doccheck-things/)

### Explanation of symbols


 Medical device

 CE marking

 Manufacturer data (address)

 Production date

 Item number

 Batch

 Unique Device Identifier

 Contains no latex/natural rubber

 Store in a dry place!


 Separate disposal of electrical and electronic appliances.

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