emergency

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These assembly instructions describe direct connection to a 12 V on-board power supply or how to fit an alternative 12 V plug to the following charging interface/connecting cable:

WM number	Charging interface/connecting cable for
WM 22436	BASE-STATION 1 NG BASE-STATION 3 NG BASE-STATION II BASE-STATION mini II
WM 10650	ACCUVAC Pro ACCUVAC Lite ACCUVAC Rescue ACCUVAC Basic
WM 28356	MEDUMAT Transport MEDUCORE Standard protective carrying case

Prerequisites for an alternative 12 V plug

- Suitable for a 13.8 V on-board power supply
- The temperature range must correspond to the internal temperatures arising during use
- Current load > 6 A
- Strain relief is essential

Remove existing 12 V plug


1. Cut connecting cable approx. 1 cm behind the 12 V plug.
2. Dispose of the plug you have cut off as electronic waste.

Connect the cable to a 12 V on-board power supply or attach the alternative 12 V plug

1. Remove the insulation from the cable ends.
2. Twist cable shielding and black cable together, then crimp with a cable end-sleeve.
3. Twist red and brown cable together, then crimp with a cable end-sleeve.
4. Connect the black cable/cable shielding to the negative contact and the red-and-brown cable to the positive contact of the 12 V on-board power supply or of the alternative 12 V plug.

Test the connection to the 12 V on-board power supply or alternative 12 V plug

1. Measure the electrical passage of the following connections:

- positive contact of connecting block/plug against positive contact of the 12 V on-board power supply or of the alternative 12 V plug
- negative contact of connecting block/plug against negative contact of the 12 V on-board power supply or of the alternative 12 V plug

2. Measure both contacts on the BASE STATION connecting block / ACCUVAC connector / MEDUMAT Transport connector and protective carrying case for MEDUCORE Standard for short-circuiting.

There may not be any short circuit.

3. If an alternative 12 V plug is used: connect the alternative 12 V plug to the on-board power supply.
4. Check the polarity of the voltage at the contacts of the BASE STATION connecting block / ACCUVAC connector / MEDUMAT Transport connector and protective carrying case for MEDUCORE Standard.

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Medical Technology GmbH + Co. KG**