
seca 336 / 336 i
seca 234 / 232 n

| | | |
|-----------|--|------------|
| Ⓓ | Gebrauchsanweisung und Garantieerklärung | 3 |
| ⒼⒷ | Instructions for use and guarantee | 39 |
| ⒻⓇ | Mode d'emploi et garantie | 75 |
| ⒾⓉ | Istruzioni per l'uso e garanzia | 111 |
| ⒺⓈ | Instrucciones de uso y garantia | 147 |
| ⓅⓉ | Manual de instruções de utilização e declaração de garantia | 183 |
| ⓅⓁ | Instrukcja obsługi i gwarancja | 219 |

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
NOTE:

This device complies with Part 15 of the FCC Rules and with RSS-210 of Industry Canada. Operation is subject to the following two conditions:

- This device may not cause harmful interference.
- This device must accept any interference received, including interference that may cause undesired operation.

NOTE:

Changes or modifications made to this equipment not expressly approved by seca may void the FCC authorization to operate this equipment.

NOTE:

Radiofrequency radiation exposure information:

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. This equipment should be installed and operated with minimum distance of 1 m between the radiator and your body. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

CONTENTS

| | | | |
|--|-----------|---|-----------|
| 1. Description of device | 40 | Taring off additional weight (Tare) | 57 |
| 1.1 Intended use | 40 | Keeping the measured result in the display (Hold) | 57 |
| Scale | 40 | Switching off the scale | 58 |
| Measuring station | 40 | 5.2 Measuring length (optional) | 58 |
| 1.2 Description of function | 40 | Measuring with the analog measuring rod ... | 59 |
| Scale | 40 | Measuring with the digital measuring rod ... | 59 |
| WiFi interface | 40 | 5.3 Measuring with barcode detection (336 i) | 59 |
| Measuring rod | 40 | Switching off the scale | 60 |
| USB interface | | 5.4 Other functions (menu) | 61 |
| (models with a WiFi interface) | 40 | Navigating in the menu | 61 |
| Data transmission | | Activating the Autohold (AHOLD) function ... | 62 |
| (models with a WiFi interface) | 40 | Activating beeps (bBEEP) | 62 |
| 1.3 User qualification | 40 | Setting filtering (FIL) | 63 |
| Administration | | Setting display backlighting (Lcd) | 63 |
| (models with a WiFi interface) | 40 | Setting standby time (AOff) | 63 |
| Operation | 40 | Calibrating the digital measuring rod (LCAL) . | 64 |
| 2. Safety information | 41 | Switching unit of length and graduations of the digital measuring rod (LUnIT) | 64 |
| 2.1 Safety precautions in these instructions for use | 41 | Restoring factory settings (rESET) | 64 |
| 2.2 Basic safety precautions | 41 | Resetting WiFi settings (rESET) (336 i) | 65 |
| Handling the device | 41 | 6. Hygiene treatment | 65 |
| Preventing electric shock | 42 | 6.1 Cleaning | 65 |
| Preventing injuries and infections | 42 | 6.2 Disinfecting | 65 |
| Preventing damage to device | 43 | Scale | 65 |
| Handling measured results | 43 | Measuring rod (if present) | 66 |
| Handling packaging material | 43 | 6.3 Sterilizing | 66 |
| Handling batteries and rechargeable batteries | 44 | 7. Function check | 66 |
| 3. Overview | 45 | 8. What do I do if | 67 |
| 3.1 Controls | 45 | 8.1 Troubleshooting | 67 |
| 3.2 Controls for seca 234 measuring rod (optional) | 46 | 8.2 Data transmission (models with a WiFi interface) | 68 |
| 3.3 Controls for seca 232 n measuring rod (optional) | 47 | 8.3 Changing batteries | 68 |
| 3.4 Menu structure (336 i) | 48 | 9. Servicing/verification | 69 |
| 3.5 Menu structure (336) | 49 | 9.1 Information on servicing and verification | 69 |
| 3.6 Symbols in the display | 50 | 9.2 Checking the content of the verification counter | 69 |
| 3.7 Markings on the device and on the type plate | 50 | 10. Technical Data | 70 |
| 3.8 Markings on the packaging | 51 | 10.1 General technical data | 70 |
| 4. Making the device operational | 52 | 10.2 Weighing technology data | 71 |
| 4.1 Scope of delivery (336 i) | 52 | 10.3 Technical data, length measurement | 71 |
| 4.2 Scope of delivery (336) | 53 | 11. Optional accessories (336) | 71 |
| 4.3 Setting up and aligning the device | 53 | 12. Optional accessories (336 i) | 72 |
| 4.4 Establishing the power supply | 54 | 13. Spare parts (336 i) | 72 |
| Inserting batteries | 54 | 14. Disposal | 72 |
| 4.5 Connecting the barcode scanner | 55 | 14.1 Disposing of the scale | 72 |
| 4.6 Setting up the WiFi connection (336 i) | 55 | 14.2 Disposing of batteries | 72 |
| Connecting the device to the WiFi network (WPS) | 55 | 15. Warranty | 72 |
| Connecting the device to the WiFi network (seca connect 103) | 56 | 16. Declaration of conformity | 73 |
| 5. Operation | 56 | 16.1 For Europe | 73 |
| 5.1 Weighing (without barcode recognition) | 56 | 16.2 For USA and Canada | 73 |
| Starting the weighing procedure | 56 | | |

1. DESCRIPTION OF DEVICE

1.1 Intended use

Scale The **seca 336/336 i** electronic baby scale is mainly used in hospitals, medical practices and inpatient care facilities in accordance with national regulations. The baby scale is for conventional determination of weight and establishment of general state of nutrition; it assists the physician supervising treatment in making a diagnosis or deciding on a course of treatment. To make an accurate diagnosis, however, the physician needs to commission other specific examinations and take their results into account, in addition to determining a weight.

Measuring station In conjunction with measuring rods **seca 234** and **seca 232 n** available as options, baby scale **seca 336/336 i** serves as a measuring station. The measuring station is for conventional determination of weight and length, establishment of general state of nutrition and assists the physician supervising treatment in making a diagnosis or deciding on a course of treatment. To make an accurate diagnosis, however, the physician needs to commission other specific examinations and take their results into account, in addition to determining weight and length.

1.2 Description of function

Scale Weight is recorded using four load cells. The scale can be transported using the carrying handle.

WiFi interface The **seca 336 i** model can transmit measured results to a hospital information system (HIS) or to seca software via a WiFi connection.

Measuring rod The scale can be expanded into a measuring station by means of analog measuring rod **seca 232 n** or digital measuring rod **seca 234**. Digital measuring rod **seca 234** has its own display and is connected to the scale with a connector. Measured results can be sent to the scale and transmitted to an HIS via WiFi.

USB interface (models with a WiFi interface) A barcode scanner can be connected to the USB interface. This allows users and patients to be identified by means of their barcodes and measured results assigned in the HIS.

Data transmission (models with a WiFi interface) The link from the device to the HIS is made via the **seca connect 103** configuration software. A current version of the configuration software can be found in the Download area for the device at www.seca.com.

1.3 User qualification

Administration (models with a WiFi interface) The device may only be set up and incorporated in a network by experienced administrators or hospital technicians.

Operation The device may only be operated by healthcare professionals.

2. SAFETY INFORMATION

2.1 Safety precautions in these instructions for use

DANGER!

Used to identify an extremely hazardous situation. If you fail to take note of this information, serious irreversible or fatal injuries will occur.

WARNING!

Used to identify an extremely hazardous situation. If you fail to take note of this information, serious irreversible or fatal injuries may result.

CAUTION!

Used to identify a hazardous situation. If you fail to take note of this information, minor to moderate injuries may result.

NOTICE!

Used to identify possible incorrect usage of the device. If you fail to take note of this information, you may damage the device, or the measured results may be incorrect.

NOTE

Includes additional information about use of the device.

2.2 Basic safety precautions

Handling the device

- ▶ Please take note of the information in these instructions for use.
- ▶ Keep the instructions for use in a safe place. The instructions for use are a component of the device and must be available at all times.

DANGER!

Risk of explosion

Do not use the device in an environment in which one of the following gases has accumulated:

- oxygen
- flammable anesthetics
- other flammable substances/air mixtures

CAUTION!

Patient hazard, damage to device

- ▶ Additional devices which are connected to electrical medical devices must provide evidence of compliance with the relevant IEC or ISO standards (e.g. IEC 60950 for data-processing devices). Furthermore, all configurations must comply with the requirements of standards for medical systems (see IEC 60601-1-1 or Section 16 of the 3rd edition of IEC 60601-1 respectively). Anyone connecting additional devices to electrical medical devices is considered a system configurer and is therefore responsible for ensuring that the system complies with the requirements of standards for systems. Your attention is drawn to the fact that local laws take precedence over the above-mentioned requirements of standards. In the event of any queries, please contact your local specialist dealer or Technical Service.
- ▶ Have servicing and verification carried out regularly as described in the relevant section of this document.
- ▶ Technical modifications may not be made to the device. The device does not contain any parts for servicing by the user. Only have servicing and repairs performed by an authorized seca Service partner. You can find service partners in your area at www.seca.com or by sending an e-mail to service@seca.com.
- ▶ Only use original seca accessories and spare parts, otherwise seca will not grant any warranty.

**CAUTION!****Patient hazard, malfunction**

- ▶ Keep other electrical medical devices, e.g. high-frequency surgical devices, a minimum distance of approx. 1 meter away to prevent incorrect measurements or wireless transmission interference.
- ▶ Keep HF devices such as cell phones a minimum distance of approx. 1 meter away to prevent incorrect measurements or wireless transmission interference.
- ▶ The actual transmission output of HF equipment may require minimum distances of more than 1 meter. Details can be found at www.seca.com.

Preventing electric shock**WARNING!****Electric shock**

- ▶ Set up the device so that the power supply socket is easy to reach and the device can be disconnected from the power supply quickly.
- ▶ Ensure that your local power supply matches the information on the power supply unit.
- ▶ Do not touch the power supply unit with wet hands.
- ▶ Do not use extension cables or power strips.
- ▶ Make sure that cables are not pinched or damaged by sharp edges.
- ▶ Make sure that cables do not come into contact with hot objects.
- ▶ Do not operate the device at an altitude of more than 3000 m above sea level.

Preventing injuries and infections**WARNING!****Injury from falling**

Baby scales generally stand on raised work surfaces. If the baby falls from this surface, this may result in serious irreversible or fatal injuries.

- ▶ Ensure that the device is steady and level.
- ▶ Route connecting cables (if present) in such a way that people cannot trip over them.
- ▶ Never leave a baby unsupervised.

**WARNING!****Risk of infection**

- ▶ Before and after every measurement, wash your hands to reduce the risk of cross-contamination and nosocomial infections.
- ▶ Hygienically reprocess the scale regularly as described in the respective section in this document.
- ▶ Make sure that the patient has no infectious diseases.
- ▶ Make sure that the patient has no open wounds or infectious skin alterations, which may come into contact with the device.

Preventing damage to device

NOTICE!

Damage to device

- ▶ Ensure that no liquids enter the device. They can damage the electronics.
- ▶ Switch off the device before disconnecting the power supply unit from the mains socket.
- ▶ For devices with power supply operation: Disconnect the power supply unit from the mains socket if you intend to not use the device for a longer period of time. Only this way it can be ensured that the device is de-energized.
- ▶ For devices with battery or rechargeable battery operation: If you do not use the device for an extended period of time, remove batteries or rechargeable batteries. Only this way it can be ensured that the device is de-energized.
- ▶ Make sure not to drop the device.
- ▶ Do not expose the device to any impacts or vibrations.
- ▶ Perform function controls regularly as described in the relevant section in this document. Do not operate the device if it is damaged or not working properly.
- ▶ Ensure that there is no heat source in the immediate vicinity. Do not expose to direct sunlight. The excessive temperature could damage the electronics.
- ▶ Avoid rapid temperature fluctuations. When the device is transported so that a temperature difference of more than 20 °C occurs, it must stay turned off for at least 2 hours before it can be turned on again. Otherwise, condensation water will form which can damage the electronics.
- ▶ Use the device only in the ambient conditions outlined in "Intended use".
- ▶ Store the device only in the storage conditions outlined in "Intended use".
- ▶ Use only chlorine and alcohol-free disinfectants which are explicitly suitable for acrylic sheet and other sensitive surfaces (active ingredient: quaternary ammonium compounds, for example).
- ▶ Do not use aggressive or abrasive cleaning agents.
- ▶ Do not use organic solvents (e.g. white spirit or petroleum spirit).

Handling measured results

NOTICE!

Inconsistent measuring results

- ▶ Before you electronically save measurement values determined using this device and use them further (e.g. in seca PC software or in a hospital information system), make sure that the measurement values are plausible.
- ▶ If measurement values are transmitted to seca PC software or a hospital information system, make sure prior to further use that the measurement values are plausible and are assigned to the correct patient.

Handling packaging material

WARNING!

Risk of suffocation

Packaging material made of plastic foil (bags) is a choking hazard.

- ▶ Keep packaging material out of reach of children.
- ▶ In the event that the original packing material may not be available anymore, only use plastic bags with security holes in order to reduce the risk of suffocation. Use recyclable materials if possible.

NOTE

Keep the original packing material for future use (e.g. returning for servicing).

Handling batteries and rechargeable batteries**WARNING!****Personal injury as a result of improper handling**

Batteries and rechargeable batteries contain harmful substances which may explode if not handled properly.


- ▶ Do not try to recharge batteries.
- ▶ Do not expose (rechargeable) batteries to heat.
- ▶ Do not burn (rechargeable) batteries.
- ▶ If acid is leaking out, avoid contact with the skin, eyes and mucous membranes. Rinse affected areas with plenty of clean water and seek medical help at once.

NOTICE!**Damage to device and malfunctions with improper handling**

- ▶ Only use the type of (rechargeable) battery specified in this document.
- ▶ When replacing (rechargeable) batteries, always replace a complete set at a time.
- ▶ Do not short-circuit (rechargeable) batteries.
- ▶ If you do not use the device for a long period of time, remove the batteries (incl. rechargeable batteries). This prevents acid from leaking into the device.
- ▶ If acid leaked into the device, discontinue use. Have the device checked by an authorised seca Service partner and repaired if necessary.

3. OVERVIEW


3.1 Controls

| No. | Control | Function |
|-----|---|--|
| 1 |  | Enter key <ul style="list-style-type: none"> • During weighing (seca 336 i only): <ul style="list-style-type: none"> - Transmit measured results to a hospital information system - Transmit measured results to seca software • In the menu: <ul style="list-style-type: none"> - Confirm selected menu item - Set value |
| 2 |  | hold arrow key <ul style="list-style-type: none"> • During weighing: <ul style="list-style-type: none"> - Activate hold function • In the menu: <ul style="list-style-type: none"> - Switch to the next menu item or select the next setting |
| 3 |  | tare arrow key <ul style="list-style-type: none"> • During weighing: <ul style="list-style-type: none"> - Activate tare function • In the menu: <ul style="list-style-type: none"> - Switch to the previous menu item or select the previous setting |
| 4 |  | Start key <ul style="list-style-type: none"> - Switch the device on and off |
| 5 | Display | Display element for measured results and to configure the device |


| No. | Control | Function |
|-----|---|--|
| 6 | Status LED | <p>Indicates the status of data recording and data transmission (seca 336 i only):</p> <ul style="list-style-type: none"> • Illuminated in green: Identification or measuring procedure started. • Flashing green (approx. 5 seconds): Measured results being sent to the HIS (depending on setting). • Illuminated in green (approx. 5 seconds): Measured results successfully sent to the HIS (depending on setting). • Illuminated in red (approx. 5 seconds, see “Data transmission (models with a WiFi interface)” from page 68): <ul style="list-style-type: none"> - ID not found in hospital information system (HIS) or in seca software. - Device has not saved measured results to the clipboard. - Measured results have not been sent to an HIS or to seca software. <p>NOTE The data recorded and transmitted are set in the configuration software. If you have any queries, please contact your administrator or hospital technician.</p> |
| 7 | Foot screw | 4 pcs, for precise alignment |
| 8 | Electrical connection (with protective cap) | For connecting a digital measuring rod |
| 9 | Battery compartment | Takes type AA batteries, 1.5 V |
| 10 | USB interface | For connecting a barcode scanner (model with WiFi interface) |
| 11 | Power supply connection | For connecting the power supply unit |
| 12 | Spirit level | Indicates whether the device is horizontal (verified models) |

3.2 Controls for seca 234 measuring rod (optional)


| No. | Control | Function |
|-----|--------------------------|--|
| 1 | Head stop | For positioning the baby |
| 2 | Connector | For connecting to a scale |
| 3 | Display | Display element for measured results |
| 4 | hold key | Sends measured result to the connected scale |
| 5 | Measuring rod with scale | For reading off length |
| 6 | Foot stop | For positioning the baby |
| 7 | Handle | For moving the foot stop |

3.3 Controls for seca 232 n measuring rod (optional)


| No. | Control | Function |
|-----|--------------------------|--------------------------|
| 1 | Head stop | For positioning the baby |
| 2 | Measuring rod with scale | For reading off length |
| 3 | Foot stop | For positioning the baby |
| 4 | Handle | For moving the foot stop |

3.4 Menu structure (336 i)


3.5 Menu structure (336)

* If digital measuring rod fitted



3.6 Symbols in the display

| | Symbol | Meaning |
|---|-------------|---|
| A | | Batteries are low |
| B | | Flashing: Select memory (BMIF) |
| C | | Differential value function (BMIF) active |
| D | | Operation with power supply unit |
| E | | Non-verifiable function active |
| F | HOLD | Hold function active |
| D | NET | Tare function active |
| E | | Select memory (BMIF function) |
| F | | Measured result in memory 1 (BMIF function) |
| G | AP | Access Point function active (service function) |
| H | | Indicates the status of the WiFi connection and the field strength of the WiFi network (models with WiFi module): <ul style="list-style-type: none"> • When setting up the WiFi network: <ul style="list-style-type: none"> - Flashing: Device connecting to a router - Illuminated: Device and router are connected and the WiFi connection has been made • In operation: <ul style="list-style-type: none"> - Illuminated: The number of bars indicates the field strength of the WiFi network |

3.7 Markings on the device and on the type plate

| Text/symbol | Meaning |
|-----------------|---|
| | Name and address of manufacturer, date of manufacture |
| REF | Model number |
| SN | Serial number |
| Mat. no. | Variant number |
| ProdID | Product identification number |
| MAC | MAC address |
| | Follow instructions for use |
| | Electrical medical device, Type B |

| Text/symbol | Meaning |
|---|--|
|  | Insulated device, protection class II |
| e | Value in mass units used to classify and verify a scale (verified models) |
| d | For electronic scales: Value which gives the difference between two consecutive display values For mechanical scales: Value in mass units giving the difference between the values of two adjacent part-marks |
|  | Class III scale to 2014/31/EU |
|  | Device complies with EC directives. <ul style="list-style-type: none"> • M: Conformity label in compliance with directive 2014/31/EU for non-automatic scales (verified models) • 16: (Example: 2016) Year in which conformity verification was performed and the CE label was applied (verified models) • 0102: Appointed office for metrology (verified models) • 0123: Appointed office for medical devices |
|  | Symbol of the US Federal Communications Commission (FCC) |
| FCC ID | Device license number from the Federal Communications Commission (FCC) |
| IC | Device license number from Industry Canada |
|  | Rating plate on the power supply connection socket <ul style="list-style-type: none"> • x-y V: Required supply voltage • max xx A: Maximum current consumption • + - : Note polarity of device plug • == : Operate device with direct current |
|  | USB interface (models with a WiFi interface) |
|  | Do not dispose of device with household waste |
|  | Device meets the requirements of the USA and Canada. Certified and tested by a licensing laboratory (NRTL) of TÜV SÜD Product Services GmbH. |


3.8 Markings on the packaging

| | |
|---|--|
|  | Protect from moisture |
|  | Arrows indicate top of product. Transport and store in an upright position. |
|  | Fragile Do not throw or drop. |
|  | Permitted min. and max. temperature for transport and storage |
|  | Permitted min. and max. moisture for transport and storage |

| | |
|---|--|
|  | Not sterile |
|  | Do not reuse |
|  | Open packaging here |
|  | Packaging material can be disposed of through recycling programs |


4. MAKING THE DEVICE OPERATIONAL

4.1 Scope of delivery (336 i)

| | Component | Pcs. |
|----------|---|------|
| a | Baby scale | 1 |
| b | Power supply unit with adapters (depending on model: power supply unit with Euro connector) | 1 |
| c | Type AA batteries, 1.5 V | 6 |
| | Instructions for use, not shown | 1 |

4.2 Scope of delivery (336)

| | Component | Pcs. |
|---|---------------------------------|------|
| a | Baby scale | 1 |
| b | Type AA batteries, 1.5 V | 6 |
| | Instructions for use, not shown | 1 |

4.3 Setting up and aligning the device

The scale is fully assembled upon delivery.

NOTICE!

Incorrect measurement as a result of force shunt

If the scale and its housing are in contact with something, e.g. a towel, weight will not be measured correctly.

► Set up the scale so that only its feet are in contact with the ground.

1. Place the scale on a firm, level surface.

NOTICE!

Incorrect alignment may cause incorrect measurements

The spirit level is very sensitive. Additional weights, such as towels, can result in poor scale alignment.


► Align the device only without a load on it.

2. Undo the knurled wheels.
3. Align the scale by turning the foot screws.
The air bubble in the spirit level must be located in the exact center of the circle.
4. Tighten the knurled wheels in the direction of the arrow.
The foot screws are secured against being adjusted.

NOTE

The scale should be checked every time it is moved to ensure it is aligned, and corrected if required.

5. If present, fit one of the measuring rods available as an option as described in the relevant assembly instructions.

4.4 Establishing the power supply


Use the batteries supplied for the following operating modes:

- Scale is used on a mobile basis
- External interfaces (e.g. WiFi, USB) are not present or are deactivated

NOTE

If you want to use the USB port or the WiFi interface (variant-dependent), the power consumption of the device rises significantly. In these cases, operate the device with the power supply unit provided.

Inserting batteries

1. Open the battery compartment.
2. Place the batteries in the battery compartment.

NOTE

Ensure the correct polarity of the batteries (markings on the battery holder). If the display **bAtt** appears, you have inserted one of the batteries the wrong way round or the batteries have discharged. If batteries are inserted the wrong way round, they must be removed again immediately.

3. Close the battery compartment.

Connecting the power supply unit

Use the power supply unit provided (variant-dependent) for the following operating modes:

- Scale is used on a stationary basis
- External interfaces (e.g. WiFi, USB) are active

NOTE

The power supply unit is available as an accessory for variants without a WiFi interface.


WARNING!

Personal injury or damage to the device as a result of incorrect power supply units


Conventional power supply units may deliver a higher voltage than is indicated on them. The scale may overheat, catch fire, melt or short-circuit.

- ▶ Use only original seca power supply units with controlled 12 V output voltage.

1. Plug the device connector required for your power supply into the power supply unit.
2. Plug the power supply connector of the power supply unit into the connecting socket of the scale.
3. Plug the power supply unit into a power supply socket.

4.5 Connecting the barcode scanner

A barcode scanner can be connected to the USB interface (variant-dependent).

Depending on the settings, the barcode scanner allows patient and user IDs to be read and the patient data called up in an HIS or in seca software. Measured results are then assigned to the patient data automatically and can be saved in the HIS.

NOTE

If you want to use the USB port or the WiFi interface (variant-dependent), the power consumption of the device rises significantly. In these cases, operate the device with the power supply unit provided.

NOTE

Observe the maximum permitted current consumption of the barcode scanner (see "General technical data" from page 70). You can find a list of recommended scanners at www.seca.com.

- ▶ Plug the USB connector of the barcode scanner into the USB socket of the device.

4.6 Setting up the WiFi connection (336 i)

Setting the WiFi interface

rF

545

1. Select the menu item "rF" (see "Navigating in the menu" on page 61).
2. Confirm your selection.
3. Select the menu item "SYS".
4. Confirm your selection.
The current setting is displayed.
5. Select the desired setting:

| Option | Effect |
|--------|---|
| On | WiFi interface automatically active in power supply mode, switches off automatically in battery mode. |
| Perm | WiFi interface remains active in both power supply and battery modes. |
| Off | WiFi interface deactivated |

6. Confirm your selection.
You exit the menu automatically.

Connecting the device to the WiFi network (WPS)

Connect your device to the WiFi network via WPS if you have access to the router.

NOTE

As soon as the device is connected to a WiFi network, the Autohold function is activated automatically. The Autohold function cannot be deactivated if the device is connected to a WiFi network.

1. Select the menu item "rF" (see "Navigating in the menu" on page 61).
2. Activate the WPS function of your router as described in the instructions for use for your WiFi router.

rF

3. Select the “Lrn” menu item.
The  symbol flashes.
The device connects to the router of the WiFi network.
As soon as the device is connected to the WiFi network, the  symbol is permanently on.

Connecting the device to the WiFi network (seca connect 103)

Connect your device to the WiFi network via the **seca connect 103** software if you do not have any access to the WPS function of the router or want to incorporate several devices:

NOTE

As soon as the device is connected to a WiFi network, the Autohold function is activated automatically. The Autohold function cannot be deactivated if the device is connected to a WiFi network.

1. Connect the barcode scanner to the device (see “Connecting the barcode scanner” on page 55).
2. Enter the data for the WiFi network in the **seca connect 103** software as described in the administrator manual for the software.
3. Scan the QR code generated in the software.
The  symbol flashes.
The name of the network and the password for the WiFi network are saved in the device.
The device connects to the router of the WiFi network.
As soon as the device is connected to the WiFi network, the  symbol is permanently on.

5. OPERATION

5.1 Weighing (without barcode recognition)

WARNING!

Injury from falling

Baby scales generally stand on raised work surfaces. If the baby falls from this surface, this may result in serious irreversible or fatal injuries.

- ▶ Ensure that the device is steady and level.
- ▶ Route connecting cables (if present) in such a way that people cannot trip over them.
- ▶ Never leave a baby unsupervised.

Starting the weighing procedure

NOTICE!

Incorrect measurement as a result of force shunt

Operation of the measuring rod affects the weight display of the scale. Weight values which are displayed during a length measurement do not correspond to the patient's actual weight.

- ▶ Ensure that you do not touch the scale whilst weight is being measured.
 - ▶ Only read off weight values before or after measuring length.
1. Ensure that there is no load on the scale.
 2. Press the Start key.

All the elements of the display are shown briefly, then **seca** appears in the display.

The scale is operational when **0.000** appears in the display.

If the scale is operated with a power supply unit, the  symbol appears in the display.

If a digital measuring rod is connected to the scale, the measuring rod is switched on automatically (see "Measuring with the digital measuring rod" from page 59).

NOTE

When the Autohold function is activated, the weight value is automatically displayed permanently (see "Activating the Autohold (AHOLD) function" on page 62).

- Place the baby on the scale.
- Press the **hold** arrow key briefly.
The weight value is displayed permanently. The  symbol (non-verifiable function) and the "HOLD" message are displayed.
- Read off the measured result.

Taring off additional weight (Tare)

Using the TARE function, you can prevent any additional weight (e.g. a towel or a cushioned mat on the weighing surface) from affecting the weighing result.

NOTICE!

Incorrect measurement as a result of force shunt

If an additional weight, e.g. a large towel, contacts the surface on which the scale is placed, weight will not be measured correctly.

► Make sure that additional weights are only placed on the scale's weighing surface.

- Place the additional weight on the scale.
- Hold down the **tare** arrow key until the "NET" message appears in the display.
- Wait until the display stops flashing and **0.000** appears instead.
- Place the baby on the scale.
- Read off the measured result.
The additional weight is deducted automatically.
- To deactivate the TARE function, press the **tare** arrow key until the "NET" message is no longer displayed or switch off the scale.

NOTE


The maximum weight which can be displayed is reduced by the weight of the objects already placed on the scale.

Keeping the measured result in the display (Hold)

When you activate the HOLD function, the weight value continues to be displayed after the load is removed from the scale. This allows you to take care of the baby before recording his or her weight.

NOTE

When the Autohold function is activated, the weight value is automatically displayed permanently (see "Activating the Autohold (AHOLD) function" on page 62).

- Place the baby on the scale.
- Press the **hold** arrow key briefly.
The display flashes until a stable weight is measured. The weight is then continuously displayed. The  symbol (non-verifiable function) and the "HOLD" message are displayed.
- To deactivate the HOLD function, press the **hold** arrow key briefly.

The  symbol and the "HOLD" message are no longer displayed. **0.000** appears in the display. If required, you can perform another weighing operation.

Switching off the scale ▶ Press the Start key.

NOTE

In battery mode, the scale remains on standby for 60 seconds. If there is no input during standby time, the scale switches off automatically. Standby time can be adjusted (see "Setting standby time (AOff)" on page 63).

5.2 Measuring length (optional)

WARNING!

Injury from falling

Baby scales generally stand on raised work surfaces. If the baby falls from this surface, this may result in serious irreversible or fatal injuries.

- ▶ Ensure that the device is steady and level.
- ▶ Route connecting cables (if present) in such a way that people cannot trip over them.
- ▶ Never leave a baby unsupervised.

CAUTION!

Injury from crushing

The baby's hands and feet can be crushed by the moving parts of the measuring rod.

- ▶ When moving the measuring rod, ensure that neither the baby's hands nor feet are resting on the measuring rod.

NOTICE!

Incorrect measurement as a result of force shunt

Operation of the measuring rod affects the weight display of the scale. Weight values which are displayed during a length measurement do not correspond to the patient's actual weight.

- ▶ Ensure that you do not touch the scale whilst weight is being measured.
- ▶ Only read off weight values before or after measuring length.


NOTICE!

Damage to device from incorrect handling

The measuring rod is not a carrying aid. The measuring rod can become deformed, difficult to move or come off the scale as a result of high strain.

- ▶ Pick up the scale only by the carrying handle if you are transporting it.

Measuring with the analog measuring rod

You can determine length using the **seca 232 n** analog measuring rod obtainable as an option.

1. Place the baby on the scale so that his or her head is in contact with the head stop of the measuring rod.
2. Pick up the measuring rod by the handle.
3. Slide the foot stop up to the baby's feet.

NOTE

Ensure that the baby's legs are stretched out straight.

4. Read off the measured result.


NOTICE!

Data loss

Analog measured values cannot be entered into the device.

- ▶ Enter analog measured values directly in the patient file to prevent data loss or the incorrect assignment of measured results.

Measuring with the digital measuring rod

You can determine length using the **seca 234** measuring rod available as an option and transmit this to the scale.

If your scale is equipped with a WiFi interface, you can send length and weight to an HIS or to seca software.

1. Place the baby on the scale so that his or her head is in contact with the head stop of the measuring rod.
2. Pick up the measuring rod by the handle of the foot stop.
3. Slide the foot stop up to the baby's feet.

NOTE

Ensure that the baby's legs are stretched out straight.

4. The measured value is shown in the display of the measuring rod.
5. Press the **hold** key on the measuring rod.
The measured value for length is transmitted to the scale.
6. In order to transmit measured values to an HIS or to seca software, proceed as described in the section entitled "Measuring with barcode detection (336 i)" from page 59.

5.3 Measuring with barcode detection (336 i)

The WiFi interface of the device connects the scale wirelessly to an HIS or to seca software.

With a barcode scanner, you can record the IDs of the baby and the user. The patient data for the baby are called up from the HIS and linked to the measured results. The measured results are then transmitted wirelessly to the HIS or to seca software.

The following conditions must be met in order to be able to identify both user and baby by means of barcodes:

- WiFi interface of the device is active
- Device is connected to an HIS
- Barcode scanner is connected to the device

NOTE

Observe the maximum permitted current consumption of the barcode scanner (see "General technical data" from page 70). You can find a list of recommended scanners at www.seca.com.

NOTE

Follow the administrator manual for the configuration software and for the barcode scanner. seca Service will be pleased to assist if you have any questions.

1. Press the Start key.
All the elements of the display are shown briefly, then **seca** appears in the display.
The scale is operational when **0.000** appears in the display.
If the scale is operated with a power supply unit, the  symbol appears in the display.
If a digital measuring rod is connected to the scale, the measuring rod is switched on automatically (see "Measuring with the digital measuring rod" from page 59).
2. Place the baby on the scale.
The baby's weight is displayed.
Wait until the status LED is illuminated green.
3. Scan your user ID.
Your user ID has been scanned successfully if "Id:U" appears in the display.
4. Scan the baby's patient ID.
The patient ID has been scanned successfully if "Id:P" appears in the display.

NOTICE!

Data loss

If there is no input for approx. 10 minutes, partial results are discarded.

▶ Measure one baby's weight and length immediately consecutively.

5. Measuring the baby's length (in conjunction with **seca 234** only, see "Measuring with the digital measuring rod" from page 59).
6. Press Enter to send measured values to the HIS or to seca software.
The status LED flashes green (approx. 5 seconds) while the data are being transmitted.
7. Wait until the status LED is illuminated green.
All measured results are present and identification is complete.
Once the measured values have been sent successfully, the status LED goes out.

NOTE

Which measured values are sent depends on the presets in the **seca connect 103** software. If you have any queries, please contact your administrator or hospital technician.

Switching off the scale

- ▶ Press the Start key.

NOTE


In battery mode, the scale remains on standby for 60 seconds. If there is no input during standby time, the scale switches off automatically. Standby time can be adjusted (see "Setting standby time (AOff)" on page 63).

5.4 Other functions (menu)

Other functions are available to you in the menu for the scale. This enables you to configure the scale to suit your conditions of use perfectly.

A summary of the menu structure can be found in “Menu structure (336 i)” on page 48.

Navigating in the menu

1. Switch on the scale.
2. Press both arrow keys simultaneously.
The last menu item selected appears in the display.
3. Press one of the arrow keys until the desired menu item appears in the display.
4. Confirm your selection by pressing Enter.
The current setting for the menu item or a submenu are displayed.
5. To amend the setting or call up a different submenu, press one of the arrow keys until the desired setting is displayed.
6. Confirm the setting by pressing Enter.
You exit the menu automatically.
7. To make more settings, call up the menu again and proceed as described.

NOTE


If no key is pressed for a few seconds, you exit the menu automatically.

Using the differential value function (bMIF)

Using the BMIF function (BMIF: Breast Milk Intake Function), you can determine the quantity of food a baby is taking during a feed. To do this, save the current weight. After the feed, call up the saved weight again and weigh the baby again. The scale calculates the difference, i.e. the quantity of food taken.

NOTE

- The weight difference is not transferred along with the WiFi data transmission.
- No measured values can be sent via WiFi if the BMIF function is activated. Ensure that the BMIF function is deactivated if you want to send measured values via WiFi (see “Data transmission (models with a WiFi interface)” on page 68).

1. Switch on the scale.
2. Place the baby on the scale before a feed.
The current weight is displayed.
3. Select the item “bMIF” from the menu.
The symbol for the BMIF function appears.
4. Confirm your selection.
5. Select the “Sto” menu item.
The \diamond symbol appears.
6. Confirm your selection.
The display “Sto 1” appears.
The \diamond symbol appears.
Arrows flash in the display.
7. Select one of three memories (in this case: 2).
8. Confirm your selection.
The current weight value is saved.
9. Place the baby back on the scale after the feed.
10. Select the item “bMIF” from the menu.

CALC

2650.0

0070.0

Activating the Autohold (AHOLd) function

AHOLd

On

Activating beeps (bEEP)

bEEP

PrESS

On

11. Confirm your selection.
12. Select the "CALC" menu item.
13. Confirm your selection.
14. Select the memory to which you saved the initial weight of the baby.
Arrows flash in the display.
15. Confirm your selection.
The weight difference, i.e. the amount of breast milk taken, is displayed.
16. To deactivate the function, select the item "bMIF" from the menu again.
17. Confirm your selection.
The function is deactivated.
You exit the menu automatically.

When you activate the Autohold function, the measured result for every weighing operation continues to be displayed after the weight has been removed from the scale. It is then no longer necessary to activate the Hold function manually for each individual weighing operation.

NOTE

- This function is activated at the factory on some models. You can deactivate the function if desired.
- As soon as the device is connected to a WiFi network, the Autohold function is activated automatically. The Autohold function cannot be deactivated if the device is connected to a WiFi network (see "Setting up the WiFi connection (336 i)" from page 55).

1. Select the item "AHOLd" from the menu.
2. Confirm your selection.
The current setting is displayed.
3. Select the desired setting:
 - On
 - Off
4. Confirm your selection.
You exit the menu automatically.

You can set whether a beep is heard every time a key is pressed and when a stable weight value is achieved. The latter is significant for the Hold/Autohold function.

NOTE

The "Beep when weight is stable" function is activated at the factory. You can deactivate this function if desired.

1. Select the item "bEEP" from the menu.
2. Confirm your selection.
3. Select a menu item:
 - PrESS: beep when key is pressed
 - HOLd: beep when weight value is stable
4. Confirm your selection.
The current setting is displayed.
5. Select the desired setting:
 - On
 - Off
6. Confirm your selection.
You exit the menu automatically.
7. If you also want to activate beeps for the second function, repeat the process.

Setting filtering (FIL)

Filtering (FIL) allows you to reduce interference (e.g. as a result of patient movement) when determining weight.

1. Select the item "FIL" from the menu.
2. Confirm your selection.
The current setting is displayed.
3. Select a filtering level.
 - 0: low filtering
 - 1: moderate filtering
 - 2: high filtering
4. Confirm your selection.
You exit the menu automatically.

Setting display backlighting (Lcd)

You can change the duration and brightness of display backlighting.

1. Select the item "Lcd" from the menu.
2. Confirm your selection.
3. Select a menu item:
 - dUr: Duration
 - brl: Brightness
4. Confirm your selection.
The current setting is displayed.
5. Select the desired setting:

| Function | Setting |
|------------|---|
| Duration | <ul style="list-style-type: none"> • Short (approx. 15 sec.) • Long (approx. 45 sec.) • PErM (permanently) |
| Brightness | <ul style="list-style-type: none"> • 50 % • 100 % • Off |

6. Confirm your selection.
You exit the menu automatically.
7. If you also want to make settings for the second function, repeat the process.

Setting standby time (AOff)

In battery mode, the scale remains on standby for 60 seconds (factory setting). If there is no input during standby time, the scale switches off automatically. Measured results will be discarded. You can adjust standby time:

| Function | Setting |
|----------------------|--------------------------|
| Factory setting | 60 seconds |
| Increment | 30 seconds |
| Minimum standby time | 30 seconds |
| Maximum standby time | 600 seconds (10 minutes) |

Proceed as follows to adjust standby time:

1. Select the item "AOff" from the menu.
2. Confirm your selection.
The time currently set is displayed (in this case: 600 seconds).
3. Select the desired time:
 - ▶ To increase value: Press **hold** arrow key
 - ▶ To reduce value: Press **tare** arrow key

4. Confirm your selection.
You exit the menu automatically.

Calibrating the digital measuring rod (LCAL)

If the value displayed in the measuring rod display deviates from the value shown on the scale of the measuring rod, adjust the measuring rod:

1. Adjust the measuring rod manually to 50 cm.
2. Compare the length to the display on the measuring rod.
3. If a different value is shown on the display, select the item "LCAL" from the menu.
4. Confirm your selection.
The calibration length of 50 cm is displayed.
5. Confirm your selection.
The measuring rod is calibrated.
You exit the menu automatically.

Switching unit of length and graduations of the digital measuring rod (LUNIT)

You can select the unit in which length is displayed in the measuring rod display. For the unit "cm", you can also adjust the graduations shown in the display.

NOTE

Observe nationally-applicable regulations governing metrology units.

1. Select the item "LUNIT" from the menu.
2. Confirm your selection.
The current setting is displayed.
3. Select the desired setting:
 - 0.1 cm
 - 0.5 cm
 - 1 cm
 - inches
4. Confirm your selection.

You exit the menu automatically.

LUNIT

CM 0.5

InCh

Restoring factory settings (rESET)

You can restore the factory settings for the functions below.

| Function | Factory setting |
|---|-----------------|
| Unit of weight | kg |
| Unit of length | cm |
| Autohold (AHOLd) | Model-dependent |
| Beep (PrESS) | Off |
| Beep (HoLd) | On |
| Filtering (FIL) | 0 |
| Standby time until the device switches off automatically (AOff) | 60 seconds |

1. Select the item "rESET" from the menu.
2. Confirm your selection.
You exit the menu automatically.
3. Switch off the scale.
Factory settings are restored and are available when the scale is switched back on.

NOTE

If you restore the factory settings, the WiFi settings will be retained. To reset the WiFi settings, proceed as described in the section entitled "Resetting WiFi settings (rESET) (336 i)" on page 65.

rESET

Resetting WiFi settings (rESET) (336 i)

You must reset the WiFi settings first before performing any reconfiguration work. The following information will be deleted in the process:

- Network name (SSID)
- Network code

NOTE

If you reset the WiFi settings, the factory settings for the scale will also be restored automatically.

1. Select the item "rESET" from the menu.
2. Press Enter until you hear a beep.
The WiFi settings are reset.
3. Set up a new WiFi connection as described in the section entitled "Connecting the device to the WiFi network (WPS)" on page 55 or "Connecting the device to the WiFi network (seca connect 103)" on page 56.

6. HYGIENE TREATMENT

WARNING! Electric shock

The device is not de-energized when the Start key is pressed and the display goes out. Use of fluids on the device may cause electric shock.

- ▶ Ensure that the device is switched off before performing any hygiene treatment.
- ▶ Disconnect the power supply connector before performing any hygiene treatment.
- ▶ Before each hygiene treatment, take the rechargeable battery out of the device (if present and removable).
- ▶ Ensure that no fluids penetrate the device.

CAUTION! Damage to device

Inappropriate detergents and disinfectants may damage the sensitive surfaces of the device.

- ▶ Use only disinfectants free of chlorine and alcohol which are explicitly suitable for acrylic sheet and other sensitive surfaces (active ingredient: quaternary ammonium compounds, for example).
- ▶ Do not use caustic or abrasive detergents.
- ▶ Do not use organic solvents (e.g. white spirit or petroleum spirit).

6.1 Cleaning

- ▶ Use a soft cloth dampened with mild soapsuds to clean the surfaces of the device.

6.2 Disinfecting

Scale

1. Check that your disinfectant is suitable for sensitive surfaces and acrylic sheet (active ingredient: quaternary ammonium compounds, for example).
2. Follow the instructions for use for the disinfectant.

3. Disinfect the device:
 - ▶ Moisten a soft cloth with disinfectant and wipe down the device.
 - ▶ Comply with the intervals: see table.

| Interval | Component |
|---------------------------------|-------------------------------------|
| Before every measurement | Tray |
| After every measurement | Tray |
| If required | Display elements, controls, housing |

Measuring rod (if present)

1. Check that your disinfectant is suitable for sensitive surfaces and acrylic sheet (active ingredient: quaternary ammonium compounds, for example).
2. Follow the instructions on the disinfectant.
3. Disinfect the device:
 - ▶ Moisten a soft cloth with disinfectant and wipe down the device with it.
 - ▶ Comply with the intervals, see table.

| Interval | Component |
|---------------------------------|---|
| Before every measurement | Head stop, foot stop |
| After every measurement | Head stop, foot stop |
| If required | <ul style="list-style-type: none"> • Controls • Housing • Measuring rod with scale |

6.3 Sterilizing

This device may not be sterilised.

7. FUNCTION CHECK

- ▶ Perform a function check prior to each use.

A complete function check includes:

- visual inspection for mechanical damage
- checking the alignment of the device
- visual and function check of the display elements
- function check of all the controls shown in the section entitled "Overview"
- function check of optional accessories

If you notice any faults or deviations during the function check, first try to resolve the error with the aid of the section entitled "What do I do if ..." in this document.

**CAUTION!****Personal injury**

If you notice any faults or deviations during the function check which cannot be resolved with the aid of the section entitled "What do I do if ..." in this document, you may not use the device.

- ▶ Have the device repaired by seca Service or by an authorized service partner.
- ▶ Follow the section entitled "Servicing/verification" in this document.

8. WHAT DO I DO IF ...


8.1 Troubleshooting

| Problem | Cause/solution |
|---|--|
| ... no weight is displayed with a load on the scale? | The device has no power supply. <ul style="list-style-type: none"> - Check whether the scale is switched on - Check whether batteries are inserted (devices operated by battery) - Check whether there is a power supply (devices operated by power supply) |
| ... 0.000 does not appear before weighing? | A load was already on the device before it was switched on. <ul style="list-style-type: none"> - Remove the load from the scale - Switch off the scale, then switch back on again |
| ... one segment is illuminated either continuously or not at all? | The corresponding point has a fault. <ul style="list-style-type: none"> - Inform seca Service |
| ... the display  appears? | Battery voltage is dropping. <ul style="list-style-type: none"> - Replace batteries |
| ... the display bAtt appears? | The batteries have discharged. <ul style="list-style-type: none"> - Replace batteries |
| ... the message StOP appears? | Maximum capacity has been exceeded. <ul style="list-style-type: none"> - Remove the load from the scale |
| ... the display tEMP appears? | The ambient temperature of the device is too high or too low. <ul style="list-style-type: none"> - Set up the device in an ambient temperature between +10 °C and +40 °C. - Wait for around 15 minutes until the device has adapted to ambient temperature. |
| ... the display Err:11 appears? | The scale has too high a load or too high a load in one corner. <ul style="list-style-type: none"> - Remove the load from the scale or distribute the load more evenly - Restart the scale - Inform seca Service |
| ... the display Err:12: appears? | The scale has been switched on with too high a load. <ul style="list-style-type: none"> - Remove the load from the scale - Restart the scale |
| ... the display Err:16 appears? | The scale was caused to oscillate and was unable to determine the zero point. <ul style="list-style-type: none"> - Restart the scale |
| ... the display Err:32 appears? | An internal communication error has occurred. <ul style="list-style-type: none"> - Restart the scale - Inform seca Service |

8.2 Data transmission (models with a WiFi interface)

| Problem | Cause/solution |
|---|---|
| ... the status LED does not come on? | The status LED is faulty. - Inform seca Service |
| ... only the "SYS" item is visible in the "rF" menu? | The WiFi interface is deactivated. - Activate WiFi interface |
| ... after calling up the menu, the "rF" item is not displayed? | The scale's WiFi interface is faulty. - Inform seca Service |
| ... the display noChG appears when the WiFi interface is activated? | Scale in battery mode, On option selected in rF\SYS menu - Select PERM option - Use power supply unit (recommended) |
| ... the status LED is illuminated red? | <ul style="list-style-type: none"> • ID not found in hospital information system (HIS) or in seca software (see fault Id:Err). • Device has not saved measured results to the clipboard. - Repeat measurement • Measured results have not been sent to an HIS or to seca software. - Deactivate BMIF function (see "Using the differential value function (bMIF)" on page 61) - Repeat measurement - Check WiFi connection |
| ... the display Id:Err appears? | The ID of the user or the patient has not been found in the HIS or in the seca software. - Scan the barcodes for the user and patient again - Inform seca Service |
| ... Enter is pressed and the display Err:71 appears? | No data transmission possible, WiFi interface is deactivated. - Activate WiFi interface (see "Setting up the WiFi connection (336 i)" on page 55) |

8.3 Changing batteries

You need 6 AA-type batteries, 1.5 V. To establish the power supply, proceed as follows:

1. Open the battery compartment.
2. Take the discharged batteries out of the battery compartment.
3. Put the new batteries in the battery compartment.

NOTE

Ensure the correct polarity of the batteries (markings on the battery holder). If the display **bAtt** appears, you have inserted one of the batteries the wrong way round or the batteries have discharged. If batteries are inserted the wrong way round, they must be removed again immediately.

4. Close the battery compartment.

9. SERVICING/VERIFICATION

9.1 Information on servicing and verification

We recommend having your device serviced prior to verification.

NOTICE!

Incorrect measurements as a result of poor servicing

- ▶ Have servicing and repairs carried out exclusively by seca Service or by an authorized service partner.
- ▶ You can find service partners in your area at www.seca.com or by sending an e-mail to service@seca.com.

Have an authorized technician perform verification according to national legal regulations.

Verification is necessary whenever one or more verification seals are damaged or the contents of the verification counter no longer match the number on the applicable verification counter sticker. If verification seals are damaged, contact seca Service directly.

9.2 Checking the content of the verification counter


This seca scale is verified. Verifications may only be performed by authorized agencies. To guarantee this, the scale is equipped with a verification counter to record any change in verification-related data.

If you want to check whether the scale has been properly verified, proceed as follows:

1. If necessary, turn off the scale.
2. Press and hold any key and start the scale.

The current content of the verification counter flashes in the display for a few seconds.

3. Compare the content of the verification counter output with the number given on the calibration counter sticker.

Both numbers have to match for the verification to be valid. If the sticker and the verification counter do not match, the scale must be verified. Please contact your service partner or seca Service. Once the scale has been verified, a new, updated verification counter sticker is used to identify the status of the verification counter. The person authorized to perform the verification secures this sticker with an additional seal. The verification counter sticker can be ordered from seca Service.

10. TECHNICAL DATA

10.1 General technical data

| General technical data | |
|--|--|
| Dimensions <ul style="list-style-type: none"> • Depth • Width • Height | 308 mm 650 mm 110 mm |
| Net weight <ul style="list-style-type: none"> • Scale without measuring rod • Scale with digital measuring rod seca 234 • Scale with analog measuring rod seca 232 n | approx. 3.7 kg approx. 5.2 kg approx. 5 kg |
| Ambient conditions, operation <ul style="list-style-type: none"> • Temperature • Air pressure • Humidity | +10 °C to +40 °C / +50 °F to 104 °F 700 - 1060 hPa 30 % - 80 %, no condensation |
| Ambient conditions, storage <ul style="list-style-type: none"> • Temperature • Air pressure • Humidity | -10 °C to +65 °C / +14 °F to 149 °F 700 - 1060 hPa 0 % - 95 %, no condensation |
| Ambient conditions, transport <ul style="list-style-type: none"> • Temperature • Air pressure • Humidity | -10 °C to +65 °C / +14 °F to 149 °F 700 - 1060 hPa 0 % - 95 %, no condensation |
| Height of digits | 20 mm |
| Power supply <ul style="list-style-type: none"> • Power supply unit (model-dependent) <ul style="list-style-type: none"> - Power supply - Maximum current consumption • Batteries <ul style="list-style-type: none"> - Power supply - Battery type | 12 V Type 500 mA 9 V 6 x type AA, 1.5 V |
| USB module (model-dependent) <ul style="list-style-type: none"> • Current consumption of scanner | Max. 500 mA |
| Supply voltage | 100 V - 240 V |
| Power supply frequency | 50 Hz - 60 Hz |
| Current consumption <ul style="list-style-type: none"> • Excluding measuring rod, excluding WiFi interface, excluding USB interface, excluding background lighting • Including digital measuring rod seca 234, WiFi interface activated, USB interface activated, permanent background lighting (brightness: 100 %) | 26 mA 240 mA |
| Maximum runtime in battery mode <ul style="list-style-type: none"> • Excluding WiFi interface, excluding USB interface, excluding background lighting • WiFi interface activated, USB interface activated | approx. 60 hours power supply unit recommended |
| Medical device in accordance with Directive 93/42/EEC | Class I with measuring function |
| EN 60601-1: <ul style="list-style-type: none"> • Insulated device, protection class II: • Electrical medical device, Type B: |   |
| Type of protection | IP20 |
| Duty cycle | Continuous duty |

17-10-07-649-100a/04-2018S

| General technical data | |
|--|---|
| Interfaces (model-dependent): | USB 2.0 |
| <ul style="list-style-type: none"> Barcode scanner seca software and hospital information system | WiFi (IEEE 802.11b/g/n/e/i) (seca connect 103 configuration software can be downloaded in the Download area for the device at www.seca.com) |

10.2 Weighing technology data

| Weighing technology data | |
|--|----------------------------|
| Verification in line with Directive 2014/31/EU | Class III |
| Maximum capacity | |
| <ul style="list-style-type: none"> Partial weighing range 1 Partial weighing range 2 | 10 kg 20 kg |
| Minimum capacity | 100 g |
| Graduations | |
| <ul style="list-style-type: none"> Partial weighing range 1, 0 kg - 10 kg Partial weighing range 2, > 20 kg | 5 g 10 g |
| Tare range | 10 kg |
| Accuracy on initial verification | |
| <ul style="list-style-type: none"> Partial weighing range 1, 0.0 kg - 2.5 kg Partial weighing range 1, 2.5 kg - 10 kg Partial weighing range 2, 10 kg - 20 kg | ± 2.5 g ± 5 g ± 10 g |

10.3 Technical data, length measurement

| Metrology data, length measurement | |
|--|--|
| Digital measuring rod seca 234 | |
| <ul style="list-style-type: none"> - Measuring range - Graduations | 35 - 80 cm 1 mm (Display graduations can be adjusted: 1 mm, 5 mm, 10 mm) |
| <ul style="list-style-type: none"> - Accuracy | ± 5 mm |
| Analog measuring rod seca 232 n | |
| <ul style="list-style-type: none"> - Measuring range - Graduations - Accuracy | 35 - 80 cm 1 mm ± 5 mm |

11. OPTIONAL ACCESSORIES (336)

| Accessory | Article number |
|---|--|
| Power supply: | |
| <ul style="list-style-type: none"> Switch-mode power supply: 100-240 V~ / 50-60 Hz / 12 V= / 0.5 A | 68 32 10 270 |
| Measuring rods: | |
| <ul style="list-style-type: none"> Analog measuring rod seca 232 n <ul style="list-style-type: none"> - cm variant - inches variant Digital measuring rod seca 234 | 232 1717 008 232 1817 008 234 1717 009 |
| Carrying bag for seca 428 baby scales | 428 0000 004 |

12. OPTIONAL ACCESSORIES (336 i)

| Accessory | Article number |
|--|--|
| Measuring rods: | |
| • Analog measuring rod seca 232 n | |
| - cm variant | 232 1717 008 |
| - inches variant | 232 1817 008 |
| • Digital measuring rod seca 234 | 234 1717 009 |
| Software (PC and server): | |
| • seca connect 103 | Download at www.seca.com |
| Barcode scanner | See recommendation at www.seca.com |
| Carrying bag for seca 428 baby scales | 428 0000 004 |

13. SPARE PARTS (336 i)

| Spare part | Article number |
|--|----------------|
| Power supply: | |
| • Switch-mode power supply: 100-240 V~ / 50-60 Hz / 12 V= / 0.5 A | 68 32 10 270 |

14. DISPOSAL

14.1 Disposing of the scale

Do not dispose of the device with household waste. The device must be disposed of properly as electronic waste. Comply with the national provisions applicable in your country. For further information contact our service department at:

service@seca.com

14.2 Disposing of batteries

Spent (rechargeable) batteries should not be discarded with household waste, regardless of whether they contain harmful substances or not. As a consumer you are obliged by law to dispose of (rechargeable) batteries via the collection points set up by the municipal authorities or the retail sector. Only discard (rechargeable) batteries when fully discharged.

15. WARRANTY

We offer a two-year warranty from the date of delivery for defects attributable to faulty material or poor workmanship. This excludes all moveable parts such as (rechargeable) batteries, cables, power supply units, etc. Defects which are covered by the warranty shall be rectified free of charge for customers on production of the sales receipt. No further claims can be accepted. The costs of shipment in both directions shall be borne by the customer where the device is not located at the customer's premises. In the event of any damage during shipment warranty claims can only be asserted where the complete original packaging was used for shipment and the scales were secured inside in the same manner as in the original packaging. You should therefore keep all packaging.

The warranty shall become null and void where the device is opened by persons not expressly authorised to do so by seca.

In the event of a warranty issue, please contact your local seca office or the dealer from whom you ordered the product.


16. DECLARATION OF CONFORMITY

16.1 For Europe

seca gmbh & co. kg hereby declares that the product meets the terms of the applicable European directives. The unabridged declaration of conformity can be found at: www.seca.com.

16.2 For USA and Canada

NOTE:

This device complies with Part 15 of the FCC Rules and with RSS-210 of Industry Canada. Operation is subject to the following two conditions:

- This device may not cause harmful interference.
- This device must accept any interference received, including interference that may cause undesired operation.

NOTE:

Changes or modifications made to this equipment not expressly approved by seca may void the FCC authorization to operate this equipment.

NOTE:

Radiofrequency radiation exposure information:
This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. This equipment should be installed and operated with minimum distance of 1 m between the radiator and your body. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

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