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Referenzen

1. Helm RE, Klausner JD, Klemperer JD, Flint LM, Huang E. Accepted but unacceptable: peripheral IV catheter failure. *J Infus Nurs*. 2015;38(3):189–203. doi: 10.1097/NAN.0000000000000100.
2. Centers for Disease Control and Prevention. *Multistate Outbreak of Burkholderia cepacia Bloodstream Infections Associated with Contaminated Prefilled Saline Flush Syringes.* 2017. Accessed March 7, 2024 at https://www.cdc.gov/hai/outbreaks/b-cepacia-saline-flush/index.html.

3. Ziegler MJ, Pellegrini DC, Safdar N. Attributable mortality of central line associated bloodstream infection: systematic review and meta-analysis. *Infection*. 2015;43(1):29–36. doi: 10.1007/s15010-014-0689-y.

4. Hu B, Tao L, Rosenthal VD, et al. Device-associated infection rates, device use, length of stay, and mortality in intensive care units of 4 Chinese hospitals: International Nosocomial Control Consortium findings. *Am J Infect Control.* 2013;41(4):301–6. doi: 10.1016/j.ajic.2012.03.037.

5. Vital signs: central line-associated bloodstream infections. Centers for Disease Control and Prevention (CDC) website. https://www.cdc.gov/vitalsigns/pdf/2011-03-vitalsigns.pdf. Published March 1, 2011. Accessed October 18, 2017.

6. Nickel B, Gorski L, Kleidon T, et al. Infusion Therapy Standards of Practice, 9th Edition. *J Infus Nurs*. 2024;47(1S Suppl 1):S1-S285. doi:10.1097/NAN.0000000000000532.

7. Bahl A, Haddad L, Hoerauf K, Mares A, Alsbrooks K. The Clinical and Economic Burdens of Infiltration and Extravasation with Peripheral Intravenous Catheters: A Contemporary Narrative Review. *Int J Nurs Health Care Res.* 2023;6:1436. https://doi.org/10.29011/2688-9501.101436.

8. Massoomi F, Arruda J. A Closer Look at Prefilled Flush Syringe Safety. *Pharmacy Practice News*. February 15, 2018. Accessed on May 19, 2023, at https://www.pharmacypracticenews.com/Clinical/Article/02-18/A-Closer-Look-at-Prefilled-Flush-Syringe-Safety/46801?sub=C63733F385E1FD4090EBC251E2F1174E9B162D4541C3E518622EA66D7AFC67F&enl=true%3Fses%3Dogst.

9. Jewett C, Nguyen D. Timeline: How Contaminated Syringes Got From Factory to Patients. June 6, 2009. Accessed May 2, 2023 at https://www.propublica.org/article/timeline-am2pat.

10. Denton A, Bodenham A, Conquest A, et al. Royal College of Nursing. Standards for Infusion Therapy, 4th Edition. https://www.rcn.org.uk/-/media/royal-college-of-nursing/documents/publications/obselete/005704.pdf?la=en. Published December 12, 2016. Accessed June 2, 2021.

11. Loveday HP, Wilson JA, Pratt RJ, et al. epic3: national evidence-based guidelines for preventing healthcare-associated infections in NHS hospitals in England. *J Hosp Infect.* 2014;86 Suppl 1:S1-S70. doi:10.1016/S0195-6701(13)60012-2.

12. Rosenthal VD, Bat-Erdene I, Gupta D, et al. International Nosocomial Infection Control Consortium (INICC) report, data summary of 45 countries for 2012-2017: Device-associated module. *Am J Infect Control.* 2020;48(4):423-432. doi:10.1016/j.ajic.2019.08.023.

13. National Patient Safety Agency. NPSA Alert - Promoting safer use of injectable medicines (NPSA 20). Published March 28, 2007. Updated August 22, 2018. Accessed October 4, 2019 at https://www.sps.nhs.uk/articles/npsa-alert-promoting-safer-use-ofinjectionmedicinesnpsa-20-2007/.

14. The Institute for Safe Medication Practices. ISMP Safe Practice Guidelines for Adult IV Push Medications. *ISMP.* Published July 14, 2015. Accessed October 20, 2020 at http://www.ismp.org/sites/default/files/attachments/2017-11/ISMP97-Guidelines-071415-3.%20FINAL.pdf.

15. The Joint Commission. National Patient Safety Goals Effective July 2020 for the Hospital Program. *The Joint Commission.* Accessed October 20, 2020 at http://www.jointcommission.org/-/media/tjc/documents/standards/national-patientsafetygoals/2020/npsg_chapter_hap_jul2020.pdf. Published March 26, 2020.

16. International Organisation for Standardisation (ISO) 13485 standard on Medical devices – Quality management systems – Requirements for regulatory purposes.

17. International Organisation for Standardisation (ISO) 14001 standard on Medical devices – Environmental best practices – Requirements for regulatory purposes.

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©2026 BD. Alle Rechte vorbehalten. BD-173648 (05/2026), angepasst von BD-157369 (04/2026)