
BD Nexiva™ Closed IV Catheter System Product configurations and specifications



BD Nexiva™
Closed IV Catheter
System – Dual Port

BD Nexiva™ Closed IV Catheter System – Dual Port


Catalog No.	Color	Gauge	Catheter Length (in)	Catheter ID (mm)	Catheter OD (mm)	Extension Tube ID (mm)	Gravity Flow Rate (mL/min)	Maximum Power Injector Flow Rate at 300 PSI (mL/sec)*	Packaging
383530	Yellow	24	0.56	0.4953- 0.5715	0.6731- 0.7493	1.22	19	N/A	20/box 80/case
383531	Yellow	24	0.75	0.4953- 0.5715	0.6731- 0.7493	1.22	18	N/A	20/box 80/case
383532	Blue	22	1.00	0.6350- 0.7112	0.8636- 0.9398	1.22	33	3.0	20/box 80/case
383536	Pink	20	1.00	0.7874- 0.8636	1.0668-1.1430	1.65	61	4.0-5.5	20/box 80/case
383537	Pink	20	1.25	0.7874- 0.8636	1.0668-1.1430	1.65	58	4.0-5.5	20/box 80/case
383538	Pink	20	1.75	0.7874- 0.8636	1.0668-1.1430	1.65	51	4.0-5.5	20/box 80/case
383539	Green	18	1.25	0.9398-1.0160	1.2700-1.3462	1.65	84	4.0-7.0	20/box 80/case
383540	Green	18	1.75	0.9398-1.0160	1.2700-1.3462	1.65	79	4.0-7.0	20/box 80/case

BD Nexiva™ Closed IV Catheter System – Single Port

Catalog No.	Color	Gauge	Catheter Length (in)	Catheter ID (mm)	Catheter OD (mm)	Extension Tube ID (mm)	Gravity Flow Rate (mL/min)	Maximum Power Injector Flow Rate at 300 PSI (mL/sec)*	Packaging
383510	Yellow	24	0.56	0.4953- 0.5715	0.6731-0.7493	1.22	19	N/A	20/box 80/case
383511	Yellow	24	0.75	0.4953- 0.5715	0.6731-0.7493	1.22	18	N/A	20/box 80/case
383512	Blue	22	1.00	0.6350- 0.7112	0.8636-0.9398	1.22	33	3.0	20/box 80/case
383516	Pink	20	1.00	0.7874- 0.8636	1.0668-1.1430	1.65	61	4.0-5.5	20/box 80/case
383517	Pink	20	1.25	0.7874- 0.8636	1.0668-1.1430	1.65	58	4.0-5.5	20/box 80/case
383518	Pink	20	1.75	0.7874- 0.8636	1.0668-1.1430	1.65	51	4.0-5.5	20/box 80/case
383519	Green	18	1.25	0.9398- 1.0160	1.2700-1.3462	1.65	84	4.0-7.0	20/box 80/case
383520	Green	18	1.75	0.9398- 1.0160	1.2700-1.3462	1.65	79	4.0-7.0	20/box 80/case

BD Nexiva™ Closed IV Catheter System – Single Port with MaxZero™ Needle-free Connector

Catalog No.	Color	Gauge	Catheter Length (in)	Catheter ID (mm)	Catheter OD (mm)	Extension Tube ID (mm)	Gravity Flow Rate (mL/min)	Maximum Power Injector Flow Rate at 300 PSI (mL/sec)*	Packaging
383550	Yellow	24	0.56	0.4953- 0.5715	0.6731-0.7493	1.22	19	N/A	20/box 80/case
383551	Yellow	24	0.75	0.4953- 0.5715	0.6731-0.7493	1.22	18	N/A	20/box 80/case
383552	Blue	22	1.00	0.6350- 0.7112	0.8636-0.9398	1.22	33	3.0	20/box 80/case
383556	Pink	20	1.00	0.7874- 0.8636	1.0668-1.1430	1.65	61	4.0-5.5	20/box 80/case
383557	Pink	20	1.25	0.7874- 0.8636	1.0668-1.1430	1.65	58	4.0-5.5	20/box 80/case
383558	Pink	20	1.75	0.7874- 0.8636	1.0668-1.1430	1.65	51	4.0-5.5	20/box 80/case
383559	Green	18	1.25	0.9398- 1.0160	1.2700-1.3462	1.65	84	4.0-7.0	20/box 80/case
383560	Green	18	1.75	0.9398- 1.0160	1.2700-1.3462	1.65	79	4.0-7.0	20/box 80/case



BD Nexiva™
Closed IV Catheter
System – Single Port



BD Nexiva™
Closed IV Catheter
System – Single Port with
MaxZero™ Needle-free
Connector

*The catheter system has been tested at the listed flow rates, however, due to variations in add-on devices, tubing, contrast media temperature, contrast media viscosity and pressure limit settings these flow rates may not be achievable.

© 2017 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company. BD-2787 (06/17)